


許博森 黃子庭



系統介紹

系統背景

攬關節是大腿骨跟骨盆腔間的關節,由一個球型(股骨頭)及杯狀物(髋臼)所組成,正常股骨頭應處於髋臼,若出現股骨頭脫離髋臼則會出現髖關節脫位情形(如嬰幼兒**臗關節發育不良)**。



投影片3

0

https://www.sundaykiss.com/%e5%85%92%e7%ab%a5%e5%ba%b7/%e9%ab%96%e9%97%9c%e7%af%80%e7%99%bc%e8%82%b2%e4%b8%8d%e8%89%af-%e ju891128@gmail.com, 2023-08-23T03:28:34.870

系統背景

攬關節發育不良為嬰幼兒最常見的髖部疾病,具體特徵像是<u>外</u> 觀可見兩腳皮膚皺褶不對稱、大腿外張時角度受限、<u>膝蓋高度不一</u> 致,但多數在新生兒時期無明顯症狀不會疼痛及不適。隨著年紀增 長就可能出現長短腳、雙腳活動度不對稱、甚至<u>跛腳</u>的情況,越晚 診斷,錯過治療的黃金時期,治療就越困難且成效較差。







投影片4

0


https://www.forcestar.com.tw/post/%E7%99%BC%E5%B1%95%E6%80%A7%E9%AB%96%E9%97%9C%E7%AF%80%E7%99%BC%E8%82%B2%E4%B8%8D%E8%89%AF%EF%BC%iu891128@gmail.com, 2023-08-23T03:30:23.075

• 系統背景

	適用階段	優	劣
徒手檢查	皆可	方便快速	操作者經驗影響大 結果較主觀
超音波檢測	0~3個月	無輻射線且相對客觀	嬰幼兒活動度影響檢測
X光檢測	4~6個月	易觀察股骨頭骨化與髋臼的 相對位置	6個月前股骨頭不清楚 較不易檢測

系統背景

α夾角	β夾角	狀況描述
>60	55	正常
55~59	55~77	髋臼不成熟 需持續追蹤
43~49	<77	發育不良 需立即治療
43~49	>77	關節鬆脫
<43	>77	半脫位
<43	>77	脫臼




https://www.yxj.org.cn/detailPage?articleId=250077 ju891128@gmail.com, 2023-08-23T07:44:31.454 0


系統目標


由於超音波檢測會受嬰幼兒在檢測時的活動度影響,因此規劃此系統,透過感測器姿態估計畫面,搭配超音波即時影像,輔助醫師能夠準確且快速地檢測到髖關節超音波影像,進而盡快做出判斷,同時縮短檢測時間,降低嬰幼兒的不適。

在醫師操作超音波檢測的同時,系統紀錄下過程中的超音波即時影像及感測器姿態估計畫面,可做為後續研究及訓練模型之資料,以利未來研究探討。



系統概括





感測器採樣

以Arduino UNO 控制兩 顆MPU6050六軸感測 器定位並測量繞軸旋轉 的變化量



數據處理

將採集到的數據進行過濾處理並運算旋轉情形(姿態估計)




顯示姿態估計

匯出姿態估計結果 並顯示兩感測器的相對狀態

硬體設定

硬體設備

Arduino UNO R3板*1 MPU6050加速度陀螺儀傳感器*2



硬體設定 Arduino

初始校正

初始化

```
#include <Wire.h>
#include <MPU6050.h>

MPU6050 mpu1;//0x68
MPU6050 mpu2;//0x69

// Pitch, Roll and Yaw values
float pitch1 = 0,roll1 = 0,yaw1 = 0;
float pitch2 = 0,roll2 = 0,yaw2 = 0;
bool shouldTransmit = false;
```

```
void setup()
 delay(500);
 Serial.begin(115200);
                                                  運動系統範圍設定
 // Initialize MPU6050 1
 while(!mpu1.begin(MPU6050_SCALE_250DPS, MPU6050_RANGE_2G,0x68))
   Serial.println("Could not find a valid MPU6050 sensor 1(0x68), check wiring!");
   delay(500);
 // Calibrate gyroscope. The calibration must be at rest.
 mpu1.calibrateGyro();
 // Set threshold sensivty.
 mpu1.setThreshold(0):
 delay(500);
                                                 運動系統範圍設定
  // Initialize MPU6050 2
 while(!mpu2.begin(MPU6050_SCALE_250DPS, MPU6050_RANGE_2G,0x69))
   Serial.println("Could not find a valid MPU6050 sensor 2(0x69), check wiring!");
   delay(500);
 // Calibrate gyroscope. The calibration must be at rest.
 mpu2.calibrateGyro();
                                                          感測器校正
 // Set threshold sensivty.
 mpu2.setThreshold(0);
 Serial.println("start computing!");
```

硬體設定 Arduino

接收數據

設定開關

```
void loop(){
  if (Serial.available() > 0) {
    char receivedChar = Serial.read();
  if (receivedChar == 's') {
      shouldTransmit = true;
      Serial.println("Start!");
  } else if (receivedChar == 'e') {
      shouldTransmit = false;
      Serial.println("END!");
  }
}
```


使用python程式, 發送啟動及結束訊 號給arduino

```
if (shouldTransmit) {
 // Read normalized values
 Vector norm1 = mpu1.readNormalizeGyro();
 Vector norm2 = mpu2.readNormalizeGyro();
 // Calculate Pitch, Roll and Yaw
 pitch1 = norm1.YAxis * timeStep://degree
 roll1 = norm1.XAxis * timeStep;
 yaw1 = norm1.ZAxis * timeStep;
 pitch2 = norm2.YAxis * timeStep;//degree
 roll2 = norm2.XAxis * timeStep;
 yaw2 = norm2.ZAxis * timeStep;
 Serial.print("(");
 Serial.print(pitch1);
 Serial.print(",");
 Serial.print(roll1);
 Serial.print(",");
 Serial.print(yaw1);;
 Serial.print(")");
 Serial.print("(");
 Serial.print(pitch2);
 Serial.print(",");
 Serial.print(roll2);
 Serial.print(",");
 Serial.println(yaw2);
 Serial.print(")");
 delay(20);
```

採樣並換算角度輸 出給電腦

數據處理與分析

處理程序



取得數據

於linux以及windows系統中,採用不同的arduino硬體偵測,用戶可根據自身系統調適

Linux

Windows

擴展卡爾曼濾波

可對當前均值和協方差的估計進行線性化的非線性濾波器通過融合可用之測量數據,能夠過濾掉有明顯誤差的測量數據

```
class ExtendedKalmanFilter:
   def init (self, initial state, initial covariance, process noise cov, measurement noise cov):
       self.state = initial state
       self.covariance = initial covariance
        self.process noise cov = process noise cov
        self.measurement noise cov = measurement noise cov
   def predict(self):
       F = np.eye(3)
       Q = self.process noise cov
       self.state = np.dot(F, self.state)
        self.covariance = np.dot(np.dot(F, self.covariance), F.T) + Q
   def update(self, measurement):
       H = np.eye(3)
       R = self.measurement noise cov
       y = measurement - np.dot(H, self.state)
       S = np.dot(np.dot(H, self.covariance), H.T) + R
       K = np.dot(np.dot(self.covariance, H.T), np.linalg.inv(S))
        self.state = self.state + np.dot(K, y)
       self.covariance = np.dot(np.eye(3) - np.dot(K, H), self.covariance)
```

0

四元數

與歐拉角同樣為表示三維空間的旋轉,為避免歐拉角萬向鎖情形而使用。

萬向鎖為整個旋轉表示系統被限制,丟失表示維度,無法繞三軸旋轉

https://kknews.cc/news/85o84v4.html ju891128@gmail.com, 2023-08-23T16:40:13.431 0

四元數

```
def update_with_new_euler(self, roll, pitch, yaw):
    new_quaternion = self.euler_to_quaternion(roll, pitch, yaw)
    self.accumulated_quaternion = self.quaternion_multiply(new_quaternion, self.accumulated_quaternion)
    self.normalize_quaternion()
```

歐拉角->四元數

```
def euler_to_quaternion(self, roll, pitch, yaw):
    roll = np.radians(roll)
    pitch = np.radians(pitch)
    yaw = np.radians(yaw)
    cr = np.cos(roll * 0.5)
    sr = np.sin(roll * 0.5)
    cp = np.cos(pitch * 0.5)
    sp = np.sin(pitch * 0.5)
    cy = np.cos(yaw * 0.5)
    sy = np.sin(yaw * 0.5)

w = cy * cp * cr + sy * sp * sr
    x = cy * cp * sr - sy * sp * cr
    y = sy * cp * sr + cy * sp * sr
    z = sy * cp * cr - cy * sp * sr
    return np.array([w, x, y, z])
```

疊加

```
def quaternion_multiply(self, q1, q2):
    w1, x1, y1, z1 = q1
    w2, x2, y2, z2 = q2

w = w1 * w2 - x1 * x2 - y1 * y2 - z1 * z2
    x = w1 * x2 + x1 * w2 + y1 * z2 - z1 * y2
    y = w1 * y2 - x1 * z2 + y1 * w2 + z1 * x2
    z = w1 * z2 + x1 * y2 - y1 * x2 + z1 * w2
    return np.array([w, x, y, z])
```

正規化

```
def normalize_quaternion(self):
   norm = np.linalg.norm(self.accumulated_quaternion)
   if norm == 0:
       return
   self.accumulated_quaternion = self.accumulated_quaternion / norm
```

相對姿態

```
# Sensor data
pitch1, roll1, yaw1 = data_thread.data_mpu1
pitch2, roll2, yaw2 = data_thread.data_mpu2

# Sensor 1
ekf1.predict()
ekf1.update(np.array([roll1, pitch1, yaw1]))
roll1, pitch1, yaw1 = ekf1.state[0:3]
handler1.update_with_new_euler(roll1, pitch1, yaw1)
sum_roll1, sum_pitch1, sum_yaw1 = handler1.quaternion_to_euler(handler1.accumulated_quaternion)

# Sensor 2
ekf2.predict()
ekf2.update(np.array([roll2, pitch2, yaw2]))
roll2, pitch2, yaw2 = ekf2.state[0:3]
handler2.update_with_new_euler(roll2, pitch2, yaw2)
sum_roll2, sum_pitch2, sum_yaw2 = handler2.quaternion_to_euler(handler2.accumulated_quaternion)
```

感測器皆經過過濾與四元數轉換

四元數->歐拉角

```
def quaternion_to_euler(self, q):
    w, x, y, z = q
# Roll
    sinr_cosp = 2 * (w * x + y * z)
    cosr_cosp = 1 - 2 * (x**2 + y**2)
    roll = np.arctan2(sinr_cosp, cosr_cosp)
# Pitch
    sinp = 2 * (w * y - z * x)
    if abs(sinp) >= 1:
        pitch = np.sign(sinp) * np.pi / 2 # use 90 degrees if out of range else:
        pitch = np.arcsin(sinp)
# Yaw
    siny_cosp = 2 * (w * z + x * y)
    cosy_cosp = 1 - 2 * (y**2 + z**2)
    yaw = np.arctan2(siny_cosp, cosy_cosp)

return roll, pitch, yaw
```

相對姿態

```
# Relative Angle
pitch ,roll ,yaw = sum_pitch2 - sum_pitch1, sum_roll2 - sum_roll1, sum_yaw2 - sum_yaw1
pitch *= 2
roll *= 2
yaw *= 2
print("sum: ",pitch ,roll ,yaw)
```

取得兩顆感測器的相對角度變化量

旋轉矩陣

```
# Rotation Matrix
rotation = euler_to_rotation_matrix(yaw,pitch,roll)
rotation_x,rotation_y,rotation_z = rotation[0:3]
```


轉換成旋轉矩陣以適用圖形顯示

```
def euler to rotation matrix(yaw, pitch, roll):
   R z = np.array([
       [np.cos(yaw), -np.sin(yaw), 0],
       [np.sin(yaw), np.cos(yaw), 0],
       [0, 0, 1]
   R y = np.array([
       [np.cos(pitch), 0, np.sin(pitch)],
       [0, 1, 0],
       [-np.sin(pitch), 0, np.cos(pitch)]
   R_x = np.array([
       [1, 0, 0],
       [0, np.cos(roll), -np.sin(roll)],
       [0, np.sin(roll), np.cos(roll)]
   # Combined rotation matrix
   R = R_z @ R_y @ R_x
   return R
```

姿態估計


座標影像

在感測器1(新生兒身上)的座標系統中 感測器2(超音波探頭)的相對旋轉 虛線:參考軸實線:即時姿態



實作影像

影像呈現



使用者操作